

(19)



Europäisches Patentamt
European Patent Office
Office européen des brevets



(11) Publication number:

0 415 695 A1

(12)

EUROPEAN PATENT APPLICATION

(21) Application number: 90309370.6

(51) Int. Cl.⁵: **G01N 33/577, G01N 33/569,
G01N 33/58, //C12N5/12**

(22) Date of filing: 28.08.90

(30) Priority: 30.08.89 US 400643

(43) Date of publication of application:
06.03.91 Bulletin 91/10

(84) Designated Contracting States:
AT BE CH DE DK ES FR GB GR IT LI NL SE

(71) Applicant: **ELI LILLY AND COMPANY**
Lilly Corporate Center
Indianapolis Indiana 46285(US)

(72) Inventor: Maciak, Ronald Steven
6350 Sunset Lane
Indianapolis, Indiana 46260(US)
Inventor: Marder, Phillip
8209 Hoover Lane
Indianapolis, Indiana 46260(US)

(74) Representative: Tapping, Kenneth George et
al
Erl Wood Manor
Windlesham Surrey, GU20 6PH(GB)

(54) **Selective cloning for high monoclonal antibody secreting hybridoma.**

(57) The invention is a method for selectively isolating cells which secrete higher levels of monoclonal antibody than did the parental culture from which they were derived. The method is based on the membrane characteristics of the cell and requires the use of a cell sorter. At present, the invention's best mode utilizes fluorescent probes in conjunction with a fluorescence activated cell sorter. The method is based on the correlation of a high level of cell surface immunoglobulin with high monoclonal antibody secretion rates. The invention provides a predictable and rapid process for selecting high monoclonal antibody producing cells.

EP 0 415 695 A1

SELECTIVE CLONING FOR HIGH MONOCLONAL ANTIBODY SECRETING HYBRIDOMAS

This invention belongs to the field of hybridoma technology and provides a predictable and rapid process for selecting and cloning cells which secrete higher than normal levels of monoclonal antibody (MAB). The process is based on the selection of those cells which present a higher than average density of immunoglobulin on their cell surfaces.

5 Monoclonal antibodies have become major constituents of the burgeoning field of biotechnology related products. It has become apparent that for the development and general distribution of MAB diagnostic and therapeutic products, large scale production techniques will be necessary.

Although MABs can be produced *in vivo* by collection of ascitic fluid, *in vitro* mass cell culture is in many ways preferable (Samoilovich et al., 1987, Hybridoma technology: new developments of practical
10 interest. *J. Immunol. Methods* 101, 153). Due to the significant production expenses incurred in these large scale cultures, obtaining maximal yields of reactive MAB is indeed important (Velez et al., 1986, Kinetics of monoclonal antibody production in low serum growth medium. *J. Immunol. Methods* 86, 45); (Reuveny et al., 1986, Comparison of cell propagation methods for their effect on monoclonal antibody yield in fermentors. *J. Immunol. Methods* 86, 53). To date, most reports in this area describe methods for
15 improving MAB yields by optimizing media formulations and by controlling environmental culture conditions (Samoilovich et al., supra). Lacking, however, are discussions of the selection and maintenance of initial variant hybridoma starter cultures capable of secreting increased MAB concentrations.

The method for producing monoclonal antibodies was first described by Kohler and Milstein in 1975. Since that time, many variations and improvements have been made on their process (for a general review,
20 see *Methods in Enzymology*, 121, 1986). This invention is not concerned with the process of constructing new hybridomas but rather focuses on the cloning and re-selection of existing hybridoma cell lines.

Hybridoma cells result from the fusion of tetraploid myeloma cells with diploid, antigen stimulated lymphocytes. The myeloma cell partner confers the ability to grow indefinitely in culture while the genetic information which codes for antibody specificity and type emanates from the lymphocyte. The resultant
25 hybridoma cell is capable of growing indefinitely and secreting a uniformly consistent and homogeneous MAB of known specificity. Due to its polyploid nature, the resultant hybridoma is also genetically unstable. Genetic instability of hybridomas is a well documented fact which can, over time, lead to profound phenotypic changes. Quite frequently, hybridoma cultures will decrease the rate at which they secrete MAB and correspondingly increase their growth rate. The combined result of those changes causes the culture
30 gradually to lose its capacity to produce and secrete MAB while consuming ever increasing amounts of substrate which results only in greater cell mass rather than MAB. Isotype switching is another phenotypic change which can occur as a hybridoma cell line is left in continuous culture. Isotype switching occurs when small sub-populations of cells within a culture change over to a different protein backbone (isotype) while maintaining the identical antigen specificity. The most detrimental phenotypic change which a
35 hybridoma culture may undergo is to lose its capacity for growth. This phenomenon is fairly common and occurs quite rapidly when cells expel or fail to replicate a chromosome which encodes a structural protein or enzyme essential for growth. In order to avoid or circumvent some of the undesired characteristics mentioned, frequent and continual re-selection or cloning of hybridoma cultures is essential so that a genetically stable, homogeneous cell line can be maintained.

40 Hybridoma cultures are cloned using any one of three basic techniques, cloning by limiting dilution, cloning over soft agarose, or cloning by fluorescence activated cell sorting. Cloning by limiting dilution is accomplished by randomly selecting a small group of cells which is then diluted into growth medium at a very low concentration. Small aliquots of the cell suspension, which theoretically contain one or less cells, are then placed into individual growth wells where they may grow into new cultures.

45 Cloning over soft agarose is accomplished by making very dilute cell suspensions in warmed growth medium containing a suitable amount of agarose. A volume of the cell suspension, which theoretically contains ten to one hundred cells, is plated into a petri dish where it is allowed to cool and solidify. As individual cells grow and divide, colonies of cells result which are held in place by the semi-solid medium. This technique is directly analogous with pour plates in bacteriology.

50 Lastly, fluorescence activated cell sorting can facilitate the cloning of a hybridoma culture. In this case, laser light is directed at individual cells as they flow through the instrument in a stream of single cells. A light scatter pattern is generated when the dense nuclear material of the cell interferes with the path of the laser beam. Thus, cells can be selected at random based on their ability to scatter laser light. After a cell has been identified, it is deflected away from the stream of cells and is directed into a growth well by means of an automated cell cloning device. The flow method affords a distinct advantage over those

previously described in that the fluorescence activated cell sorter is capable of depositing one and only one cell into a particular growth well. The other methods rely on theoretical calculations which statistically make the chances of one cell per well or one cell per unit area likely rather than assured.

All three of the aforementioned techniques rely on the random selection of cells. The successful isolation of a clone with desirable MAb secretion characteristics is no more predictable than the laws of random probability will allow. This invention provides a method of selecting hybridoma cells for single cell cloning which are most likely to secrete at a higher level than those cloned by random selection

The invention comprises a method of selectively isolating a population of cells which secrete monoclonal antibody of known specificity at a level higher than the parent culture from which they were derived and which have immunoglobulin associated with, and presented on the outer surface of, the plasma membrane, wherein cells of the parent culture are contacted with a probe which preferentially binds to the immunoglobulin and those cells having high probe intensity are isolated by means of a cell sorter. The invention further comprises a method of individually depositing the selected cells into separate growth wells by means of a cell sorter outfitted with an automated cloning accessory. Lastly, the invention also comprises a method of rescuing the monoclonal antibody secreting capacity of an essentially non-secreting hybridoma culture.

While the following detailed discussion of the invention focuses on fluorescent probes and fluorescence activated cell sorters, other methods of detection of immunoglobulin on the cell surface are entirely usable. The following descriptions that employ fluorescent probes and fluorescence activated cell sorters are merely used to illustrate the present invention. Fluorescence detection is preferred, however.

A hybridoma cell is a eukaryotic cell which is formed by the fusion of one or more myeloma cells with one or more lymphocytes. Numerous different types of myeloma cells are commercially available through the American Type Culture Collection (Rockville, Maryland) and Ventrex Laboratories, Inc. (Portland, Maine). Lymphocytes can be obtained from a variety of organ sites located within an immunized donor. Most typically, the sites are the spleen, the lymph nodes and the blood. Examples of established hybridoma cell lines that are stored and catalogued by the American Type Culture Collection (ATCC) on which this invention could be utilized follow:

Antigenic determinant	Type	Designation	ATCC No.
Acetylcholinesterase, human	IgG ₁	AE-2	HB73
Actin, prokaryotic	IgG ₁	ACT I	HB 80

	Antigenic determinant	Type	Designation	ATCC No.
5	Bacterial adhesins (K99 pili)	IgG	2BD4E4 K99	HB 8178
	Carcinoembryonic antigen (CEA), high mw	IgM	1116NS-3d	CRL 8019
10	Choriocarcinomas	IgG ₁	162-46.2	HB 187
	Clathrin, bovine brain	IgG _{2a}	CVC.7	TIB 138
	Colorectal carcinoma	IgG ₁	1116-NS-19-9	HB 8059
	Clq, human	IgG	4A4B11	HB8327
15	Clq, human	IgG	12A5B7	HB 8328
	Diphtheria toxin	IgG (human)	16M3F10	HB 8363
20	ds DNA	IgM	CH26-1352	HB 8329
	DNA polymerase α	IgG ₁	STK 1	CRL 1652
	<u>Escherichia coli</u> (K99 pilus)	IgG	2BD4E4 K99	HB 8178
	Fibronectin	IgG ₁	P ₃ NP/PFN	HB 91
25	Forssman antigen	IgM (rat)	M1/22.25.8.HL	TIB 121
	Glutamic acid decarboxylase	IgG ₁	GAD-1	HB 184
30	HLA antigens			
	A2	IgG _{2b}	BB7.2	HB 82
	A11 and A24	IgG ₃	A11.1M	HB 164
35	A, B, C, monomorphic	IgG ₁	BB7.5	HB 120
	B7 and B4	IgG ₁	BB7.6	HB 115
	DR, DP, and DQ	IgG ₁	IVA12	HB 145
	Ia-human, monomorphic	IgG _{2a}	L243	HB 55
40	IgE, human	IgG _{2a}	E5BB3IIA2	HB 121
	IgG Fc region, human	IgG _{2a}	HP 6017	CRL 1753
	IgG _{2a} & IgG _{2b} , mouse	IgG _{2a}	ED1-19-1-6-5	HB 90
45	IgG _{2a} & IgG _{2b} , mouse	IgG _{2a} (rat)	ED1-19-1-6-5	HB90

50

55

	Antigenic determinant	Type	Designation	ATCC N .
	IgG ₃ hinge regi n, human	IgG ₁	HP 6047	CRL 1774
5	IgG ₄ Fc region, human	IgG ₁	HP 6025	CRL 1775
	IgM, Fab portion	IgG ₁	M-2E6	HB 138
	Insulin	IgG ₁	CE9 H9	HB 127
10	Insulin, human, residue A8-10	IgG ₁	AE9D6	HB 125
	y Interferon, human	IgG ₁	IFGCP-F1BA10	HB 8291
	Interleukin-4	IgG ₁	11B11	HB 188
15	<u>Legionella pneumophila</u>	IgG _{2b}	Lpl MAB 2	CRL 1770
	κ Light chains, human	IgG _{2a}	HP 6053	CRL 1758
	Light harvesting complex (LHC-II)	IgG ₁	MLH1	CRL 1766
20	Low density lipoprotein (LDL) receptors, rabbit	IgG ₁	9D9	CRL 1703
25	Lung cancer, non-small cell	IgG _{2a} κ	704A1	HB 8302
	Lung cancer, small cell	IgM	SM1	HB 8462
	Lymphocytes, activated, human	IgG ₁	OKT9	CRL 8021
30	Lyt 2.2	IgG _{2b} (rat)	2.43	TIB 210
	Mammary carcinoma	IgG ₁	B6.2	HB 8106
	Melanoma	IgG _{2a}	XMMME-001	HB 8759
35	Monoamine oxidase B, human	IgG	MAO-1C2#81	HB 8176
	c-myc protein, human	IgG ₁	MYC CT9-B7.3	CRL 1725
	v-myb protein, viral	IgG ₁	MYB 2-7.77	CRL 1724
40	Neuroblastoma cells	IgM	Mab 126	HB 8568
	Null cells, human	IgG _{2b}	OKM1	CRL 8026
	Ornithine decarboxylase,	IgM	B11	HB 8372

45

50

55

	Antigenic determinant	Type	Designation	ATCC No.
	mouse kidney			
5	Renin, hog	IgG ₁	F32 VIII C4	CRL 1653
	Ricin A chain	IgG ₁	TFTA1	CRL 1771
	Ricin B chain	IgG ₁	TFTB1	CRL 1759
10	Small cell carcinoma, human	IgM	SM1	HB 8462
	T cell subset and macrophages, swine	IgMk	76-6-7	HB 141
15	T cell surface antigen, human	IgM	2T8-3E10	HB 8213
	T cells, peripheral, human	IgG _{2a}	OKT3	CRL 8001
	Tetanus toxoid	IgG (human)	9F12	HB 8177
20	Thrombospondin	IgG _{2a}	ahTSP-1	HB 8432
	L-Thyroxine	IgG ₁	T4 clone 3	HB 8499
	Transferrin receptor	IgG _{2a}	L5.1	HB 84
25	2, 4, 6 Trinitrophenyl	IgG ₁	1B7.11	TIB 191
	Viral Antigens			
	Epstein-Barr virus	IgG ₁	72A1	HB 168
30	Hepatitis B surface antigen	IgG ₁	H25B10	CRL 8017
	Herpes simplex type 1	IgG _{2a}	39-S	HB 8180
	Influenza virus	IgG ₃	HK-PEG-1	CL 189
35	Influenza A nucleoprotein	IgG ₁	46/4	HB 67

Those skilled in the art will realize and understand that the hybridoma cell lines listed above are merely illustrative of cell lines that can be obtained or constructed according to known methods and that such list in no way limits the present invention.

Hybridoma cells produce and secrete MABs into the surrounding medium. Hybridomas also have MABs of the same type and specificity associated, either integrally or passively, with the outer surface of the plasma membrane. At present, the exact mechanism by which these MABs appear on the cell surface is unknown and in any event, the present invention does not reside in such mechanism. For the purposes of this patent, MABs which appear on the outer cell surface and which are either integrally or passively associated with the membrane will be referred to as immunoglobulin.

The present invention utilizes the correlation between high immunoglobulin expression and high MAB secretion. This correlation is high and thus by selecting hybridomas with a high density of immunoglobulin, the general likelihood of isolating a high MAB secreting culture is greatly increased over conventional methods.

The nature of the probe used to identify and select the cells of interest may vary. The probe should accomplish two functions to enable the operation of this invention. It should first preferentially bind to the immunoglobulin located on the outer surface of the cell. Secondly, it should be able to fluoresce in a manner compatible with fluorescence activated cell sorter analysis. So long as the probe is able to accomplish these two functions, its exact composition is not critical.

Those skilled in the art will appreciate that this invention is not limited to a one step attachment of the probe to the immunoglobulin. Two part or even multipartite probes are within the scope of this invention and are within the skill of the art.

For example, biotin - avidin is a well known two part probe system. In this case, an antibody or protein

which specifically binds the immunoglobulin is chemically conjugated to biotin. Such biotinylated antibodies and proteins are commercially available (see table below). Avidin, which specifically binds to biotin, is the second part of the probe and is commercially available in forms conjugated to fluorescein and Texas Red® (Boehringer Mannheim, Indianapolis, IN). The biotinylated protein, which is specific for the immunoglobulin, constitutes the first part of the probe. The second part of the probe is fluorochrome labeled avidin. Thus the fluorochrome is ultimately bound to the immunoglobulin through the specificity of the biotinylated antibody or protein.

Techniques whereby a fluorescent signal is amplified, by binding antibodies to other antibodies which have already specifically bound a target, are also well known in the field of immunology and are considered multipartite probes. An example of such a probe is constructed from goat antiserum, that is specific for mouse immunoglobulins, and is bound by sheep antiserum that is specific for goat antibodies. Then, fluorochrome tagged rabbit antibodies, that specifically bind sheep antibodies, are allowed to bind to the sheep antibodies. In this instance, the goat antiserum is the part of the probe which specifically binds the immunoglobulin located on the surface of a mouse hybridoma. The fluorochrome labeled rabbit serum is the part of the probe which fluoresces. The goat and sheep antiserum are the components of the probe which link the fluorescent portion to the probe part which specifically binds to the cell surface immunoglobulin.

A generalized probe which binds to the Fc region of any antibody produced by a given species may also be used. These types of probes are standard immunological reagents which include MAbs and antisera raised against various types of antibody classes in different host species. MAbs and antisera are purified and immunoadsorbed to insure proper specificity. The purified reagent is then chemically coupled to a standard fluorochrome such as fluorescein, rhodamine, Texas Red®, phycoerythrin or phycocyanin. Many commercial sources for this probe type exist. Some examples follow:

<u>Source</u>	<u>Specificity</u>	<u>Fluorochrome</u>	<u>Company</u>
Rat (MAb)	Mouse Kappa Chain	Fluorescein	Amac Inc.: Westbrook, ME
Mouse (MAb)	Rabbit IgG	(Biotin)	Amac Inc.: Westbrook, ME
Rabbit	Human IgG	Rhodamine	Calbiochem: San Diego, CA
Goat	Human IgG	Rodamine	Calbiochem: San Diego, CA

EP 0 415 695 A1

	<u>Source</u>	<u>Specificity</u>	<u>Fluorochrome</u>	<u>Company</u>
	Goat	Mouse IgG	Rhodamine	Calbiochem: San Diego, CA
5	Goat	Mouse IgM	Rhodamine	Calbiochem: San Diego, CA
	Goat	Rabbit IgG	Rhodamine	Calbiochem: San Diego, CA
10	Goat	Rat IgG	Rhodamine	Calbiochem: San Diego, CA
	Goat	Mouse IgG	Texas Red®	Calbiochem: San Diego, CA
15	Goat	Rat IgG	Texas Red®	Calbiochem: San Diego, CA
	Goat	Rabbit IgG	Phycoerythrin	Calbiochem: San Diego, CA
20	Goat	Human IgG	Phycoerythrin	Calbiochem: San Diego, CA
25	Goat	Mouse IgG	Phycocyanin	Calbiochem: San Diego, CA
30	Goat	Mouse IgA	(Biotin)	Calbiochem: San Diego, CA
	Sheep	Human IgG-Fc	Fluorescein	Cappel: Malvern, PA
35	Goat	Mouse IgG	Rhodamine	Cappel: Malvern, PA
	Rabbit	Human IgG	Texas Red®	Cappel: Malvern, PA
	Goat F(ab') ₂	Mouse IgG	Fluorescein	Cappel: Malvern, PA
	Goat	Human IgE	Fluorescein	Cappel: Malvern, PA
40	Rabbit	Human IgG	Fluorescein	Cappel: Malvern, PA
	Sheep	Human IgG	Fluorescein	Cappel: Malvern, PA
45				
50				
55				

	<u>Source</u>	<u>Specificity</u>	<u>Fluor chrome</u>	<u>Company</u>
	Rat (MAB)	Rabbit IgG	Rhodamine	Jansen, Inc.: Piscataway, NY
5	Mouse (MAB)	Rat IgD	Rhodamine	Jansen, Inc.: Piscataway, NY
	Rat (MAB)	Mouse IgG ₂ [®]	Fluorescein	Pandex Division: Mundelein, IL
10	Rat (MAB)	Mouse IgM	Fluorescein	Pandex Division: Mundelein, IL
15	Mouse (MAB)	Mouse Kappa Chain	Fluorescein	Pandex Division: Mundelein, IL
	Rat (MAB)	Rabbit IgG	Fluorescein	Sera Lab: Westbury, NY
20	Mouse (MAB)	Rat IgG _{2b}	Fluorescein	Sera Lab: Westbury, NY
25	Mouse (MAB)	Rat IgG _{2a}	Fluorescein	Sera Lab: Westbury, NY
	Rat (MAB)	Mouse IgG ₁	Fluorescein	Zymed: San Fran. CA
	Mouse (MAB)	Human IgG ₁	(Biotin)	Zymed: San Fran. CA
30	Mouse (MAB)	Human IgG _{2a}	(Biotin)	Zymed: San Fran. CA
	Rat (MAB)	Mouse IgG _{2a}	(Biotin)	Zymed: San Fran. CA
	Rat (MAB)	Mouse IgG _{2b}	(Biotin)	Zymed: San Fran. CA
	Rat (MAB)	Mouse IgM	(Biotin)	Zymed: San Fran. CA
35	Goat	Human IgG	(Biotin)	Zymed: San Fran. CA

Those skilled in the art will realize and understand that the probes listed above are merely illustrative of probes that can be obtained or constructed according to known methods and that such list in no way limits the present invention.

The particular antigen to which the Fab region of the immunoglobulin binds can also be used as a probe providing the antigen has been tagged with a fluorochrome or has inherent fluorescent properties. For example, if a hybridoma culture producing MABs specific for insulin were to be selectively cloned using this invention, fluorescein labeled insulin (Sigma Chemical Co., St. Louis, Mo.) could be used as the specific probe. Two other examples of antigen type probes are fluorescein labeled bovine serum albumin and fluorescein labeled dog serum albumin, both of which are commercially available (Sigma Chemical Co., St. Louis, MO.). These reagents could be used as probes when using the invention on hybridoma cultures that produce MABs specific for bovine serum albumin and dog serum albumin respectively.

Since virtually all purified proteins can be readily fluoresceinated using well established techniques, antigen type probes can be constructed which are suitable for use with the invention. Commercially available, purified, illustrative proteins such as carcinoembryonic antigen (CEA: Biodesign Inc., Kennebunkport, ME), human chorionic gonadotropin (hCG: Cambridge Medical Diagnostics Inc., Bellerica, ME), and human growth hormone (hGH: Hybritech, San Diego, CA) can be fluorochrome labeled. Once tagged with a fluorochrome, these proteins are suitable probes when utilizing the invention on hybridoma cultures which produce MABs directed against CEA, hCG and hGH respectively. Using established techniques for labeling proteins, a great variety of antigen probes can be constructed which are useful in the invention.

A hapten is a chemical compound which is unable to elicit an immunological response on its own.

However, when the compound is conjugated to a true immunogen, the immune system is able to respond and form antibodies to the previously non-immunogenic compound. One such immunogen is the keyhole limpet hemocyanin (KLH) protein. The antibiotic actinomycin-D (Act-D; Sigma Chemical Co.; St. Louis, Mo.) is an example of a haptenic substance. When the immune system is challenged with Act-D it is unable to respond. Once coupled to KLH, Act-D acts as an antigenic epitope on the highly immunogenic KLH protein. Hybridomas which secrete MABs directed towards Act-D have been constructed. Skilled artisans will understand that a fluorochrome can be coupled directly to Act-D or to a carrier protein such as KLH by conventional means and that the resultant conjugate can be used as a probe in this invention. Other hapten probes can be constructed in a similar fashion. Haptenic compounds such as 2,4-dinitrophenol and 2,4,6-trinitrophenol when coupled to fluorochrome labeled carrier molecules such as KLH are suitable hapten probes for the invention.

Many other possibilities exist for constructing hapten probes. Compounds such as phosphocholine and phosphotyrosine (Boehringer Mannheim, Indianapolis, IN) have been used as haptens to generate MABs. These compounds are suitable hapten probes for the invention when conjugated to a fluorochrome labeled carrier protein. In general, any compound capable of eliciting an immune response as a hapten can either be directly labeled with a fluorochrome or be conjugated to a fluorochrome labeled carrier protein and used as a hapten probe in the invention.

A particular antigen or hapten could possess the appropriate fluorescent characteristics to enable the use of the native antigen or hapten itself as the probe in the invention. Examples of this specialized situation would include all of the standard fluorochromes fluorescein, rhodamine, Texas Red®, phycoerythrin and phycocyanin when used as haptens.

The descriptions used to illustrate the invention use a fluorescence activated cell sorter (for a brief overview see Herzenberg et al., 1976, Fluorescence activated cell sorting. Sci. Amer. 234, 108; for an in depth analysis see Flow Cytometry and Sorting, eds. Melamad, Mullaney and Mendelsohn. John Wiley and Sons, Inc., New York, 1979). Briefly, fluorescence activated cell sorters take a suspension of cells and pass them single file by a laser placed near a detector. The detector measures the fluorescent intensity of each cell as it passes through the instrument and generates a histogram plot of cell number versus fluorescent intensity. Gates or limits can be placed on the histogram thus identifying a particular population of cells. The instrument is then instructed to electrostatically separate the selected cells away from all the other cells. In the batch mode, the fluorescence activated cell sorter will deflect those cells which fall within the gated area into a common vessel. The resultant batch collected culture is not a monoclonally derived culture and is, therefore, genetically heterogeneous at the time the selection process is completed. The heterogeneous condition does not negate the fact that the newly batch sorted culture will secrete MAB at a higher rate than did the culture from which it was derived.

Currently, fluorescence activated cell sorters can be equipped with automated cell cloning devices. Such a device enables one to instruct the instrument to singly deposit a selected cell into an individual growth well. The preferred method of operating the invention is to establish a new culture by depositing a selected single cell into a growth well where it is allowed to grow into a monoclonally derived hybridoma culture. Thus, genetic homogeneity is established within the newly cloned culture.

In an isolation according to the present invention, the cell sorter is set up to select those cells having the highest probe intensity, usually a small fraction of the cells in the culture. The level of intensity at which the sorter is set and the fraction of cells which is selected, depend on the condition of the parent culture and the objective of the isolation. In general, the operator should first sort an aliquot of the culture, and record the histogram of intensity versus number of cells. The operator can then set the selection level and isolate an appropriate number of the most active cells.

The invention is particularly useful when a culture's MAB secreting capacity has deteriorated to such a point that MAB is no longer detectable in the spent medium. The cell line is now considered an essentially non-secreting culture. This phenomenon is well known to those skilled in the art and occurs relatively often for reasons discussed in the Background section. Even though MAB levels are below the detection limits of several very sensitive assays, the culture is not totally non-secreting. In an essentially non-secreting culture, cells which retain the ability to secrete MAB may be in the range of one in 1,000 to one in 10,000. A viable, secreting culture can be isolated from an essentially non-secreting culture using conventional cloning techniques but only through very aggressive, extensive and time consuming efforts. The invention circumvents the laborious and time consuming process of randomly cloning cells which have already lost their capacity to secrete MAB. The invention enables one to sort through hundreds of thousands of cells in a matter of minutes and identify and capture those vanishingly few cells which are still able to secrete MAB.

Regardless of the cloning technique used, the resultant hybridoma cultures must be screened for MAB production and specificity. A wide variety of techniques and methods for measuring these parameters is

described in the scientific literature.

It should be stressed that the present invention is not limited to conventional hybridoma cells. Methods utilizing any cell type that produces and secretes a MAb of known specificity are within the scope of the present invention. The invention is not limited by the type of MAb produced. Methods utilizing cells which produce chimeric antibodies, MAbs of any species, bifunctional MAbs and multispecific MAbs are all within the scope of the present invention.

Example 1

L2-KS

The murine hybridoma cell line L2-KS has been described in detail (Starling et al., 1989, In vivo efficacy of monoclonal antibody-drug conjugates of three different subisotypes which bind the human tumor-associated antigen defined by the KS 1/4 antibody. *C.I.I.* 28, 171) and is available through the American Type Culture Collection under HB 9940. Briefly, it is a murine hybridoma which secretes an IgG₁ MAb that reacts with a 40,000 dalton human tumor associated glycoprotein known as KS 1/4 (Varki et al., 1984, Antigens associated with a human lung adenocarcinoma defined by monoclonal antibodies. *Cancer Res.* 44, 681.)

All flow cytometric cell sorting was done with an Epics C (Coulter Epics Division, Hialeah FL) flow cytometer-cell sorter (FCCS). Hybridoma cells were excited with 600 MW of 488 nm light emitted from an argon ion laser. The resultant cell-associated fluorescence light was captured through 488 interference and 515 long pass absorbance optical filters. The optical alignment of the FCCS Instrument was optimized and the fluorescence readings obtained were standardized using uniform fluorospheres. The flow cytometer was outfitted with an automated cloning accessory called an "AUTO-CLONE®" (Coulter Epics Division, Hialeah FL) which can be instructively programmed to dispense a single particle or cell into selected compartment of a standard 96 well microtiter culture plate. The sort delay timing was assured by microscopic observation of sorted fluorospheres. In a test evaluation of the AUTO-CLONE®, a single fluorosphere was programmed to be deposited into wells of three separate microtiter plates. Upon fluorescence microscopic inspection, each of the 288 wells had exactly 1 fluorosphere, with no doublets and no misses. The cell-contact components of the FCCS were sterilized prior to cloning by flushing with 95% ethanol for 30 min followed by exhaustive rinsing in sterile phosphate buffer saline (PBS; GIBCO; Grand Island, N.Y.) sheath fluid.

Cells were prepared for sorting by harvesting $1-3 \times 10^6$ hybridoma cells and incubating them in 1 ml of a 1:10 dilution of F(ab')₂ fragments of goat-anti-mouse-IgG-FITC (AMIGG-FITC; Cappel, Cat #1311-0081, West Chester, PA) for 45 min on ice followed by two washes in culture medium (HL-1; Ventrex Labs; Portland, MA; supplemented with 10% fetal bovine serum; Hyclone; Ogden, UT.). 200 μ l of a 50 μ g/ml solution of propidium iodide (PI, Sigma Chemical Co., St. Louis MO) was added to 1 ml of the cell preparation just prior to sorting to help define and select viable cells.

Subpopulations of cells to be sorted were FCCS defined by either rectilinear gates or by polygonal bitmaps of collected two parameter (FALS vs. log green fluorescence) histograms. A three droplet sorting procedure was employed to help assure maximum purity (per manufacturer's instructions). Single hybridoma cells were deposited into microtiter wells which were seeded 24 hr. previously with 5×10^6 C3H/HeJ thymocyte feeder cells (Jackson Labs; Bar Harbor, MA) in 200 μ L of culture medium.

MAb concentrations in the spent culture media were determined by automated laser nephelometry using a Behring Laser Nephelometer (Behring Diagnostics, La Jolla, CA.). Concentrations are expressed μ g/ml/72 hours and have been normalized to one million cells. Results appear in Table 1 which further exemplify the present invention.

Example 2

L1-KS

The murine hybridoma cell line L1-KS has been described in detail (Starling et al. supra). It secretes an IgG_{2b} MAb that reacts with the KS 1/4 antigen. Cloning was performed in substantial accordance with the teachings of Example 1. Results appear in Table 1 which further exemplify the present invention.

Example 3

L4-KS

The murine hybridoma cell line L4-KS has been described in detail (Starling et al. supra). It secretes an IgG_{2a} MAb that reacts with the KS 1/4 antigen. Cloning was performed in substantial accordance with the teachings of Example 1. Results appear in Table 1 which further exemplify the present invention.

Example 4

14-95-55

The murine hybridoma cell line 14-95-55 produces an IgG_{2a} MAb that reacts with CEA. This line was an essentially non-secreting culture at the time it was used in the invention. Cloning was performed in substantial accordance with the teachings of Example 1. Results appear in Table 1 which further exemplify the present invention.

Table 1

	Secretion Rate	Secretion Rate
	Pre-invention	Pre-invention
Hybridoma Name	ug/ml/72 hrs/1x10 ⁶ cells	ug/ml/72 hrs/1x10 ⁶ cells
L1-KS	10.9	34.1
L2-KS	8.2	18.5
L4-KS	6.8	25.1
14-95-55	<1.5	40.8

Claims

1. A method of selectively isolating a population of cells which secrete monoclonal antibody of known specificity at a level higher than the parent culture from which they were derived and which have immunoglobulin associated with, and presented on the outer surface of, the plasma membrane, wherein cells of the parent culture are contacted with a probe which preferentially binds to the immunoglobulin and those cells having high probe intensity are isolated by means of a cell sorter.
2. The method of Claim 1 wherein the culture is a hybridoma culture, the probe is fluorescent and the cell sorter is a fluorescence activated cell sorter.
3. The method of Claim 2 wherein the probe is chosen from the group consisting of fluorochrome-tagged antibody, fluorochrome-tagged antigenic substance or fluorochrome-tagged haptenic substance.
4. The method of Claim 1 wherein the selected population of cells are individually deposited into separate growth wells by means of a cell sorter outfitted with an automated cloning accessory.

5. The method of Claim 4 wherein the culture is a hybridoma culture, the probe is fluorescent and the cell sorter is a fluorescence activated cell sorter.
6. The method of Claim 5 wherein the probe is chosen from the group consisting of fluorochrome-tagged antibody, fluorochrome-tagged antigenic substance or fluorochrome-tagged haptenic substance.
- 5 7. The methods of Claims 1, 2, 3, 4, 5, and 6 wherein the parent culture is essentially non-secreting.

10

15

20

25

30

35

40

45

50

55



European
Patent Office

EUROPEAN SEARCH REPORT

Application Number

EP 90 30 9370

DOCUMENTS CONSIDERED TO BE RELEVANT			
Category	Citation of document with indication, where appropriate, of relevant passages	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int. Cl.5)
X,P	EP-A-0 347 230 (BIOGENEX LABS.) * Column 5, lines 8-26 *	1-3,6,7	G 01 N 33/577 G 01 N 33/569 G 01 N 33/58 // C 12 N 5/12
Y	NL-A-8 301 244 (NEDERLANDSE CENTRALE ORGANISATIE VOOR TOEGEPAST-NATUURWETENSCHAPPELIJK ONDERZOEK TE 'S-GRAVENHAGE) * Page 2, line 35 - page 3, line 26; page 4, line 27 - page 6, line 15 *	1-3,6,7	
Y	WO-A-8 701 385 (U. VON ARNIM) * Abstract; claim 5 *	1-3,6,7	
A	DE-A-3 806 558 (HITACHI LTD) * Column 2, line 44 - column 3, line 56 *	1-3,6	
A	TRENDS IN BIOTECHNOLOGY, vol. 6, no. 9, September 1988, pages 209-213; G.A. NEIL et al.: "Recent improvements in the production of antibody-secreting hybridoma cells" * Page 211, columns 2-3 *	1-3	
			TECHNICAL FIELDS SEARCHED (Int. Cl.5)
			G 01 N
The present search report has been drawn up for all claims			
Place of search The Hague		Date of completion of search 26 November 90	Examiner HITCHEN C.E.
<div>CATEGORY OF CITED DOCUMENTS</div> <div>X: particularly relevant if taken alone Y: particularly relevant if combined with another document of the same category A: technological background O: non-written disclosure P: intermediate document T: theory or principle underlying the invention</div> <div>E: earlier patent document, but published on, or after the filing date D: document cited in the application L: document cited for other reasons &: member of the same patent family, corresponding document</div>			